



Human Parvovirus B19 IgG antibody ELISA Kit

Enzyme Immunoassay for the qualitative and semi-quantitative determination of IgG-class antibodies to Parvovirus B19 in serum.

Catalog number: ARG80596

For research use only. Not for use in diagnostic procedures.

TABLE OF CONTENTS

SECTION	Page
PRINCIPLE OF THE ASSAY	3
MATERIALS PROVIDED & STORAGE INFORMATION	3
MATERIALS REQUIRED BUT NOT PROVIDED	4
TECHNICAL HINTS AND PRECAUTIONS	4
SAMPLE COLLECTION & STORAGE INFORMATION	5
REAGENT PREPARATION.....	5
ASSAY PROCEDURE	6
CALCULATION OF RESULTS	7
INTERPRETATION OF RESULTS	7

MANUFACTURED BY:

Arigo Biolaboratories Corporation

Address: No. 22, Ln. 227, Gongyuan Rd., Hsinchu City 300, Taiwan

Phone: +886 (3) 562 1738

Fax: +886 (3) 561 3008

Email: info@arigobio.com

Human Parvovirus B19 IgG antibody ELISA Kit ARG80596

PRINCIPLE OF THE ASSAY

This assay employs the qualitative enzyme immunoassay technique. Recombinant Parvovirus B19 antigen (VP1 proteins) has been pre-coated onto a microtiter plate. Standards or samples are pipetted into the wells and any Ab present is bound by the immobilized antigen. After washing away any unbound substances, a HRP-conjugated anti human antibody is added to each well and incubate. After washing away any unbound antibody-enzyme reagent, a substrate solution (TMB) is added to the wells and color develops in proportion to the amount of Ab bound in the initial step. The color development is stopped by the addition of acid and the intensity of the color is measured at a wavelength of 450 nm \pm 2 nm. The concentration of Ab in the sample is then determined by comparing the O.D of samples to the standard curve.

MATERIALS PROVIDED & STORAGE INFORMATION

Store the unopened kit at 2-8 °C. Use the kit before expiration date.

Component	Quantity	Storage information
Antigen-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
Negative Control	2 ml (Ready-to-use)	4°C.
Cut-Off standard	2 ml (Ready-to-use)	4°C.
Positive Control	1 ml (Ready-to-use)	4°C.
HRP-Antibody	20 ml	4°C

Human Parvovirus B19 IgG antibody ELISA Kit ARG80596

conjugate	(Ready-to-use)	
20X Wash buffer	30 ml	4°C
Sample Diluent	100 ml	4°C
TMB substrate	14 ml	4°C (Protect from light)
STOP solution	14 ml	4°C

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader capable of measuring absorbance at 450nm
- Pipettes and pipette tips
- Deionized or distilled water
- Automated microplate washer (optional)

TECHNICAL HINTS AND PRECAUTIONS

- Wear protective gloves, clothing, eye, and face protection especially while handling blood or body fluid samples.
- Store the kit at 4°C at all times.
- Briefly spin down the antibody conjugate concentrate and HRP-Streptavidin concentrate before use.
- If crystals are observed in the 20X Wash buffer, warm to RT (not more than 50°C) until the crystals are completely dissolved.
- Ensure complete reconstitution and dilution of reagents prior to use.
- It is highly recommended that the standards, samples and controls be assayed in duplicates.
- Change pipette tips between the addition of different reagent or samples.

SAMPLE COLLECTION & STORAGE INFORMATION

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Serum - Use a serum separator tube (SST) and allow samples to clot for 30 minutes before centrifugation for 15 minutes at 1000 x g. Remove serum and assay immediately or aliquot and store samples at ≤ -20 °C. Avoid repeated freeze-thaw cycles.

REAGENT PREPARATION

- **1X Wash buffer:** Dilute 20X Wash buffer into distilled water to yield 1X Wash buffer. The diluted Wash Solution is stable for 4 weeks at 2 °C to 8 °C. If crystals in the solution disappear, warming up the solution to 37 °C in a water bath. Be sure that the crystals are completely dissolved before use.
- **Patient sample:** Prior to assaying dilute each patient specimen 1+100 with Sample Diluent; e.g. 10 μ L of specimen + 1 mL of Sample Diluent, mix well, let stand for 15 minutes, mix well again before use..
Please note: Controls are ready for use and must not be diluted!

ASSAY PROCEDURE

All materials should be equilibrated to room temperature (RT) before use. Standards, samples and controls should be assayed in duplicates.

1. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal it.
2. Leave 1 well for substrate blank, 1 well for the negative control, 2 wells for the Cut-off control and 1 well for the positive control.
3. Add 100 μ l of diluted samples and controls into respective wells.
4. Incubate for 60 minutes at 37°C.
5. Aspirate each well and wash, repeating the process 4 times for a total 5 washes. Wash by filling each well with 1 \times Wash Buffer (350 μ l) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating, decanting or blotting against clean paper towels.
6. Add 100 μ l 1X Antibody solution into each well (expect the substrate blank). Incubate for 30 minutes at RT.
7. Wash as according to step 5.
8. Add 100 μ l of TMB Reagent to each well. Incubate for 15 minutes at room temperature.
9. Add 100 μ l of Stop Solution to each well.
10. Read the OD with a microplate reader at 450/629 nm immediately. It is recommended read the absorbance within 30 min after adding stop solution.

CALCULATION OF RESULTS

- Adjust the ELISA microplate or microstrip reader to zero using the substrate blank in well A1. If - due to technical reasons - the ELISA reader cannot be adjusted to zero using the substrate blank in well A1, subtract the absorbance value of well A1 from all other absorbance values measured in order to obtain reliable results!
- Measure the absorbance of all wells at 450 nm and record the absorbance values for each control and patient sample in the distribution and identification plan.
- Dual wavelength reading using 620 nm as reference wavelength is recommended.
- Where applicable calculate the mean absorbance values of all duplicates.

INTERPRETATION OF RESULTS

CO = Mean of OD of Cut-off value

POSITIVE: Patient (mean) absorbance values more than 20 % above CO
(Mean OD patient > 1.2 x CO)

GREY ZONE: Patient (mean) absorbance values from 20 % above to 10 % below CO repeat test 2 - 4 weeks later with new patient samples (0.9 x CO ≤ Mean OD patient ≤ 1.2 x CO)

Results in the second test again in the grey zone ⇒ **NEGATIVE**

Human Parvovirus B19 IgG antibody ELISA Kit ARG80596

NEGATIVE: Patient (mean) absorbance values more than 10 % below CO
(Mean OD patient < 0.9 x CO)

$$\frac{\text{Patient (mean) absorbance value} \times 10}{CO} = \text{Unit}; U$$

Cut-off value: 10 U

Grey zone: 9 - 12 U

Negative: < 9 U

Positive: > 12 U

Validation of the Test Run

The test run may be considered valid provided the following criteria are met:

Neg. Control: Absorbance value lower than 0.200

Cut-off Control: Absorbance value between 0.350 – 0.850

Pos. Control: Absorbance value between 0.650 – 3.000

Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is 100%.

Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is 98%.