



Fumonisin (Rapid) ELISA Kit

Enzyme Immunoassay for the rapid quantitative determination of
Fumonisin in cereals and beer /gyle

Catalog number: ARG80800

For research use only. Not for use in diagnostic procedures.

TABLE OF CONTENTS

SECTION	Page
INTRODUCTION	3
PRINCIPLE OF THE ASSAY	4
MATERIALS PROVIDED & STORAGE INFORMATION	4
MATERIALS REQUIRED BUT NOT PROVIDED	5
TECHNICAL HINTS AND PRECAUTIONS	5
SAMPLE COLLECTION & STORAGE INFORMATION	6
REAGENT PREPARATION.....	6
ASSAY PROCEDURE	7
CALCULATION OF RESULTS	8
EXAMPLE OF TYPICAL STANDARD CURVE	8
QUALITY ASSURANCE.....	9

MANUFACTURED BY:

Arigo Biolaboratories Corporation

Address: No. 22, Ln. 227, Gongyuan Rd., Hsinchu City 300, Taiwan

Phone: +886 (3) 562 1738

Fax: +886 (3) 561 3008

Email: info@arigobio.com

INTRODUCTION

Fumonisin in addition to zearalenone, deoxynivalenol and other trichothecenes belongs to the fusarium toxins. These toxins are already produced on the field in consequence of a contact of the cereals by fusarium species. These toxins show an extreme stability against high temperatures (up to 100°C), and they can remain active in contaminated food for years. Fumonisin can be found in maize, oats and other types of grain. Worldwide a contamination in maize of 60% has been detected. When ingested by animals, fumonisin leads to neurotoxicity, hepatotoxicity and lung edema, mainly in horses and pigs. Therapeutic measures are the change of the grain given to the animals, or the administration of diuretic drugs. In human patients hints for the appearance of esophagus cancer could be associated with the exposition to fumonisin B1. Values assessed for the acute toxicity are 8 mg per kg weight and for the chronic situation 25 mg/kg in feed stuff.

Since June 2010 the US Food and Drug Association recommends maximum amounts of 2- 100 ppm for raw cereals depending on the intended use. In the European Union the limits are 0.1 – 0.5 ppm for food and 5-60 ppm for feed products. Thus a monitoring of food and feed with respect to the concentration of fumonisin is obligatory.

The Fumonisin RAPID ELISA represents a highly sensitive detection system and is particularly capable of the rapid quantification of fumonisin contaminations in cereals and beer.

Fumonisin (Rapid) ELISA Kit ARG80800

PRINCIPLE OF THE ASSAY

This assay employs the quantitative competitive enzyme immunoassay technique. An antibody binding protein is coated on the surface of a microtiter plate. Fumonisin containing samples or standards, a fumonisin-peroxidase conjugate and an antibody directed against fumonisin are given into the wells of the microtiter plate. The conjugate competes with the fumonisin of Samples / standards for the limited number of antibody sites. Simultaneously the anti-fumonisin antibody is bound to the antibody-binding protein coated on the microtiter plate. After 10 minutes incubation at room temperature the wells are washed with diluted washing solution to remove unbound material. A substrate solution is added and incubated for 10 minutes, resulting in the development of a blue color. The color development is inhibited by the addition of a stop solution, and the color turns yellow. The yellow color is measured at 450 nm. The concentration of fumonisin is indirectly proportional to the color intensity of the test sample.

MATERIALS PROVIDED & STORAGE INFORMATION

Store the unopened kit at 2-8 °C. Use the kit before expiration date.

Component	Quantity	Storage information
Antibody-coated microplate	12 strips x 8-well	4°C
Anti-Fumonisin Antibody	6 ml (ready to use)	4°C
HRP-Fumonisin	6 ml (ready to use)	4°C
Standards (0, 0.05, 0.15, 0.5, 1.5, 5 ppm)	6 X 1 ml (ready to use)	4°C

Fumonisin (Rapid) ELISA Kit ARG80800

Sample Diluent	2 X 60 ml (ready to use)	4°C
10x Wash Buffer	60 ml	4°C
TMB substrate	15 ml	4°C (Protect from light)
STOP solution	15 ml	4°C

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader capable of measuring absorbance at 450nm
- Pipettes and pipette tips
- Deionized or distilled water
- Automated microplate washer (optional)

TECHNICAL HINTS AND PRECAUTIONS

- Wear protective gloves, clothing, eye, and face protection especially while handling blood or body fluid samples.
- Store the kit at 4°C at all times.
- Briefly spin down the anti-Fumonisin antibody and HRP-Fumonisin conjugate before use.
- If crystals are observed in the 10X Wash buffer, warm to RT (not more than 50°C) until the crystals are completely dissolved.
- Ensure complete reconstitution and dilution of reagents prior to use.
- It is highly recommended that the standards, samples and controls be assayed in duplicates.
- Change pipette tips between the addition of different reagent or samples.
- Samples contain azide cannot be assayed.

SAMPLE COLLECTION & STORAGE INFORMATION

Cereals

- Grind sample to pass through a 20 mesh sieve and thoroughly mix prior to sub-sampling.
- Suspend 20 g of sample in 100 mL of double 70 % methanol.
- Mix suspension for 5 minutes.
- Filter through Whatman #1 filter or alternatively centrifuge at a minimum of 3000 g for 5 minutes.
- Dilute 50 μ L of filtrate/supernatant with 700 μ L of sample diluent and test the sample in the ELISA.

Beer / Gyle

- Dilute an adequate volume of sample diluent with 5 % methanol.
- Carbonized beer samples should be preliminarily degassed by moderate heating.
- Cloudy beers (such as beer brewed from wheat) / gyle should preliminarily be sterile-filtered.
- Dilute 100 μ L beer / gyle with 900 μ L sample diluents/methanol dilution.

In case of too high concentrated samples, an adequate volume of sample diluent is diluted with 5 % methanol. The sample extracts have to be further diluted with this dilution.

REAGENT PREPARATION

- **1X Wash buffer:** Dilute 10X wash buffer into distilled water to yield 1X wash buffer.

ASSAY PROCEDURE

All materials should be equilibrated to room temperature (RT) before use. Standards, samples and controls should be assayed in duplicates.

1. Remove excess microplate strips from the plate frame, return them to the foil pouch containing the desiccant pack, and reseal it.
2. Add 100 μ l of standards and samples in duplicate into wells. Immediately add 50 μ l HRP-Fumonisin conjugate and 50 μ l anti-fumonisin antibody into each well.
3. Incubate for 10 minutes at RT.
4. Aspirate each well and wash, repeating the process 2 times for a total 3 washes. Wash by filling each well with 1X wash buffer (350 μ l) using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating, decanting or blotting against clean paper towels.
5. Add 100 μ l of TMB mixture to each well. Incubate for 10 minutes at room temperature in dark.
6. Add 100 μ l of Stop Solution to each well.
7. Read the OD with a microplate reader at 450 nm immediately.


Fumonisin (Rapid) ELISA Kit ARG80800

CALCULATION OF RESULTS

1. Calculate the average absorbance values for each set of standards, controls and patient samples.
2. Using linear graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred method. Other data reduction functions may give slightly different results. The diluted samples must be further converted by the appropriate dilution factor according to the sample preparation procedure as described above.

EXAMPLE OF TYPICAL STANDARD CURVE

The following data is for demonstration only and cannot be used in place of data generations at the time of assay.



QUALITY ASSURANCE

Sensitivity

The limit of detection (LOD) of the Fumonisin Rapid test is 0.015 ppm.

Validation experiments with common matrices resulted in the following LODs [ppm].

Wheat	0.010
Rye	0.025
Barley	0.011
Oats	0.015
Corn	0.014
Rice	0.025
Beer	0.004

The limit of quantification (LOQ) of the Fumonisin Rapid test is 0.05 ppm.

Due to the variety of sample matrices and their influence on the blank, results less than the LOQ should be treated as negative.

Intra-assay and Inter-assay Precision

The CV value of intra-assay precision was 3-4% and the CV of inter-assay precision was 6-8%

Recovery

Wheat flour	92%
Oats flour	89%
Corn flour	94%
Rice flour	105%
Beer flour	105%