



## **T-2 Toxin ELISA Kit**

Enzyme Immunoassay for the quantification of T-2 Toxin in food (extraction, dilution).

Catalog number: ARG81034

---

For research use only. Not for use in diagnostic procedures.

## **TABLE OF CONTENTS**

<b>SECTION</b>	<b>Page</b>
INTRODUCTION.....	3
PRINCIPLE OF THE ASSAY.....	4
MATERIALS PROVIDED & STORAGE INFORMATION.....	5
MATERIALS REQUIRED BUT NOT PROVIDED.....	6
TECHNICAL NOTES AND PRECAUTIONS.....	7
SAMPLE COLLECTION & STORAGE INFORMATION.....	8
REAGENT PREPARATION.....	9
ASSAY PROCEDURE.....	9
EXAMPLE OF TYPICAL STANDARD CURVE.....	10
CALCULATION OF RESULTS.....	11
QUALITY ASSURANCE.....	12

### **MANUFACTURED BY:**

Arigo Biolaboratories Corporation

Address: No. 22, Ln. 227, Gongyuan Rd., Hsinchu City 300, Taiwan

Phone: +886 (3) 562 1738

Fax: +886 (3) 561 3008

Email: [info@arigobio.com](mailto:info@arigobio.com)

### INTRODUCTION

T-2 Mycotoxin (pronounced as 'Tee-Two') is a trichothecene mycotoxin. It is a naturally occurring mold byproduct of *Fusarium spp.* fungus which is toxic to humans and animals. The clinical condition it causes is alimentary toxic aleukia and a host of symptoms related to organs as diverse as the skin, airway, and stomach. Ingestion may come from consumption of moldy whole grains. T-2 can be absorbed through human skin. Although no significant systemic effects are expected after dermal contact in normal agricultural or residential environments, local skin effects can't be excluded. Hence, skin contact with T-2 should be limited. [Provide by Wikipedia: T-2 mycotoxin]

### **PRINCIPLE OF THE ASSAY**

This assay employs the competitive enzyme immunoassay technique. The samples or Standards (T-2 Toxin Standards) are first added to the Antibody-binding Protein coated microplate. Then T-2 Toxin conjugate and Antibody specific for T-2 Toxin is added to each well and incubate. The T-2 Toxin Conjugate competes with the samples / Standards for the limited number of Antibody sites. After washing away any unbound substances, the TMB Substrate is added to the wells and color develops in inversely proportion to the amount of T-2 Toxin content bound with the Antibody. The color development is stopped by the addition of STOP Solution and the intensity of the color is measured at a wavelength of 450 nm. The concentration of T-2 Toxin in the sample is then determined by comparing the O.D of samples to the standard curve.

## T-2 Toxin ELISA Kit ARG81034

---

### MATERIALS PROVIDED & STORAGE INFORMATION

Upon receipt, Store the unopened kit at 2-8 °C. Use the kit before expiration date.

Component	Quantity	Storage information
Antibody-Binding Protein coated microplate	8 X 12 strips	4°C
Standards (0, 17.5, 87.5, 350, 875, 1750 ppb T-2 Toxin Standards, dyed red)	1 mL each (ready to use)	4°C
Diluent Buffer (dyed red)	60 mL X 2 (ready to use)	4°C
T-2 Toxin Conjugate (T-2 Toxin-Peroxidase, dyed red)	6 mL (ready to use)	4°C
Antibody (Anti-T-2 Toxin Antibody, dyed blue)	6 mL (ready to use)	4°C
10X Wash Buffer	60 mL	4°C
TMB Substrate	15 mL (ready to use)	4°C (protect from light)
STOP Solution	15 mL (ready to use)	4°C

### **MATERIALS REQUIRED BUT NOT PROVIDED**

- Microplate reader capable of reading at 450 nm (620 nm as optional reference wave length)
- Deionized or distilled water
- Methanol
- Vortex mixer or Ultra-Turrax
- Plastic bag to store unused microtiter strips.
- Pipettes and pipette tips
- Multichannel micropipette reservoir
- Microtiter plate washer (recommended)

### TECHNICAL NOTES AND PRECAUTIONS

- Wear protective gloves, clothing, eye, and face protection especially while handling blood or body fluid samples.
- Store the kit at 2-8°C at all times.
- Briefly spin down the all vials before use.
- If crystals are observed in the 10X Wash Buffer, warm to 37°C until the crystals are completely dissolved.
- Minimize lag time between wash steps to ensure the plate does not become completely dry during the assay.
- Ensure complete reconstitution and dilution of reagents prior to use.
- Take care not to contaminate the TMB Substrate. Do not expose the TMB solution to glass, foil or metal. If the solution is blue before use, DO NOT USE IT.
- Do NOT return leftover TMB Substrate to bottle. Do NOT contaminate the unused TMB Substrate. If the solution is blue before use, DO NOT USE IT.
- Change pipette tips between the addition of different reagent or samples.
- Plate readers measure vertically. Do not touch the bottom of the wells.
- Use a new adhesive plate cover for each incubation step.
- Taping the well strips together with lab tape can be done as an extra precaution to avoid plate strips from coming loose during the procedure.
- Include a standard curve each time the assay is performed.
- Run both standards and samples in at least duplicates (triplicate is recommended).

## **SAMPLE COLLECTION & STORAGE INFORMATION**

### **Cereals / Meat:**

- A. Grind sample to pass through a 20 mesh sieve and thoroughly mix prior to sub-sampling.
- B. Suspend 20 g of sample in 100 mL of 70% methanol.
- C. Mix suspension for 5 minutes.
- D. Centrifuge at a minimum of 3000 g for 5 minutes.
- E. Dilute 100  $\mu$ L of supernatant with 600  $\mu$ L of Diluent Buffer and test the samples in the ELISA.

### **Beer / Milk:**

- A. Depending on the number of samples. Dilute an adequate volume of Diluent Buffer with 10% methanol.
- B. Carbonized samples should be preliminarily degassed by moderate heating.
- C. Cloudy samples (such as beer brewed from wheat) should preliminarily be sterile-filtered.
- D. Dilute 100  $\mu$ L sample with 900  $\mu$ L Diluent Buffer/methanol dilution and test the diluted samples in the ELISA.

### **Note:**

1. Aliquot samples for testing and store at  $-80^{\circ}\text{C}$ . Avoid repeated freeze-thaw cycles. In case of too highly concentrated samples, an adequate volume of Diluent Buffer is diluted with methanol to a concentration of 10% methanol. The sample extracts have to be further diluted with this dilution.
2. Samples containing sodium azide should not be used in the assay.



### REAGENT PREPARATION

- **1X Wash Buffer:** Dilute 10X Wash Buffer into distilled water to yield 1X Wash Buffer. (E.g., add 50 ml of 10X Wash Buffer into 450 ml of distilled water to a final volume of 500 ml) The 1X Wash Buffer is stable for up to 4 weeks at 2-8°C.

### ASSAY PROCEDURE

All materials should be equilibrated to room temperature (RT, 20-25°C) before use. Each Samples or Standards should be assayed in duplicate or triplicate.

1. Add **100 µL** of **prepared samples or Standards** to the Antibody-Binding Protein coated microplate.
2. Add **50 µL** of the **T-2 Toxin Conjugate** to each well.
3. Add **50 µL** of the **Antibody** to each well.
4. Incubate at **RT** for **10 mins** on a microplate shaker.
5. Aspirate each well and wash, repeating the process 2 times for a total 3 washes. Wash by filling each well with **1× Wash Buffer (300 µl)** using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating, decanting or blotting against clean paper towels.
6. Warm **TMB Substrate** to **RT**. Add **100 µl** of **TMB Substrate** to each well, including the blank wells. Incubate for **10 minutes** at **RT** in the dark.
7. Add **100 µl** of **Stop Solution** to each well, including the blank wells. The color of the solution should change from blue to yellow.

## T-2 Toxin ELISA Kit ARG81034


---

- Read the OD with a microplate reader at **450 nm** immediately. (optional: read at 620 nm as reference wavelength) It is recommended reading the absorbance **within 30 minutes** after adding the stop solution.

### EXAMPLE OF TYPICAL STANDARD CURVE

The following figures demonstrate typical results with the T-2 Toxin ELISA Kit. One should use the data below for reference only. This data should not be used to interpret actual results.

T-2 Toxin (ppb)	% binding of 0 ppb
0	100
17.5	85
87.5	62
350	35
875	20
1750	13



### CALCULATION OF RESULTS

1. Calculate the average absorbance values for each set of standards and samples.
2. Using log-log, semi-log or linear graph paper, construct a standard curve by plotting the mean absorbance obtained from each standard against its concentration with absorbance value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Use the mean absorbance value for each sample determine the corresponding concentration from the standard curve.
4. Due to a deviating sample preparation process the results for Beer / Milk samples additionally have to be **multiplied with 0.286** in order to get the real concentration of the sample.
5. Automated method: The results in the IFU have been calculated automatically using a 4 PL (4 Parameter Logistics) curve fit. 4 Parameter Logistics is the preferred method. Other data reduction functions may give slightly different results.
6. arigo provides GainData®, an in-house development ELISA data calculator, for ELISA data result analysis. Please refer our GainData® website for the detail. (<https://www.arigobio.com/elisa-analysis>)
7. If the samples have been diluted, the concentration read from the standard curve must be further converted by the appropriate dilution factor according to the sample preparation procedure as described above.

## T-2 Toxin ELISA Kit ARG81034

---

### QUALITY ASSURANCE

#### Sensitivity

The limit of detection (LOD) of the T-2 Toxin ELISA kit is 13 ppb.

#### Intra-assay and Inter-assay precision

The CV value of intra-assay was 3-4 % and inter-assay was 3-6 %.

#### Recovery

Wheat	100 %
Rye	103 %
Barley	96 %
Oats	97 %
Rice	95 %
Corn	92 %
Meat	97 %
Beer	105%
Milk	95%

#### Cross-reactivity

Cross-reactivity	Relative to T-2 Toxin (=100%)
HT-2 Toxin	3.0 %
T-2 Triol	0.35 %
T-2 Tetraol	0.07 %