

Human HBV Core antigen / HBcAg IgM antibody ELISA Kit

Human HBV Core antigen / HBcAg IgM antibody ELISA Kit is an Enzyme Immunoassay kit for the determination of Human HBcAg IgM antibody in serum and plasma (EDTA, heparin, citrate).

Catalog number: ARG82307

Package: 96 wells

For research use only. Not for use in diagnostic procedures.

TABLE OF CONTENTS

SECTION	Page
INTRODUCTION	3
PRINCIPLE OF THE ASSAY	3
MATERIALS PROVIDED & STORAGE INFORMATION	4
MATERIALS REQUIRED BUT NOT PROVIDED	5
TECHNICAL NOTES AND PRECAUTIONS	5
SAMPLE COLLECTION & STORAGE INFORMATION	7
REAGENT PREPARATION	8
ASSAY PROCEDURE	8
CALCULATION OF RESULTS	9
INTERPRETATION OF RESULT	11
QUALITY ASSURANCE	11

MANUFACTURED BY:

Arigo Biolaboratories Corporation

Address: No. 22, Ln. 227, Gongyuan Rd., Hsinchu City 300, Taiwan

Phone: +886 (3) 562 1738

Fax: +886 (3) 561 3008

Email: info@arigobio.com

INTRODUCTION

HBcAg (core antigen) is a hepatitis B viral protein. It is an indicator of active viral replication; this means the person infected with Hepatitis B can likely transmit the virus on to another person.

HBcAg is an antigen that can be found on the surface of the nucleocapsid core (the inner most layer of the hepatitis B virus). While both HBcAg and HBeAg are made from the same open reading frame, HBcAg is not secreted. HBcAg is considered "particulate" and it does not circulate in the blood but recent study show it can be detected in serum by Radioimmunoassay. However, it is readily detected in hepatocytes after biopsy. The presence of both HBcAg and HBeAg proteins together act as a marker of viral replication, and antibodies to these antigens are a marker of declining replication. [Provide by Wikipedia: HBcAg]

PRINCIPLE OF THE ASSAY

This assay employs the enzyme immunoassay technique. The solid phase of the microplate is made of polystyrene wells coated with Anti-Human IgM. When a serum or plasma specimen containing Anti-HBc IgM is added to the Anti-Human IgM coated plate and incubated, IgM antibodies present in the specimen bind to the Anti-Human IgM on the wells. After addition of an HBcAg-containing reagent and a solution containing peroxidase-conjugated anti-HBc a further incubation takes place, during which (Anti-Human IgM) \cdot (HBcAg) \cdot (Anti-HBc Peroxidase) complex is formed on the wells. After washing the plate to remove unbound material, a TMB Substrate is added to the wells and

incubated. If Anti-HBc IgM is present in the specimen, after washing, the activity of peroxidase on the wells reflects the content of anti-HBc IgM in a specimen. Then the Stop Solution is added to the wells to stop the reaction. The optical density of developed color is read with a plate reader at O.D. 450 nm.

MATERIALS PROVIDED & STORAGE INFORMATION

Store all other components at 2-8°C. Do not freeze. Use the kit before expiration date.

Component	Quantity	Storage information
Anti-IgM Coated Microplate	96 wells	4°C
Anti-HBc·Peroxidase Solution	8 mL (ready to use)	4°C
Anti-HBc IgM Positive Control	2.5 mL (ready to use)	4°C
Anti-HBc IgM Negative Control	2.5 mL (ready to use)	4°C
Specimen Diluent Buffer	36 mL (ready to use)	4°C
HBcAg Reagent	8 mL (ready to use)	4°C
20X Wash Buffer	58 mL	4°C
TMB Substrate A	12 mL	4°C (protect from light)
TMB Substrate B	12 mL	4°C (protect from light)
Stop Solution	12 mL (ready to use)	4°C

Note:

- Adhesive Slips / Absorbent Pads / Black Cover provide as needed.
- All component except Anti-IgM Coated Microplate are stable for up 1 month once open. (Anti-IgM Coated Microplate is stable for up 2 month once open)

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader capable of reading at 450 nm (620-690 nm as optional reference wave length)
- Centrifuge and centrifuge tube
- Incubator (37°C)
- Deionized or distilled water
- Pipettes and pipette tips
- Multichannel micropipette reservoir
- Microtiter plate washer (recommended)

TECHNICAL NOTES AND PRECAUTIONS

- Wear protective gloves, clothing, eye, and face protection especially while handling blood or body fluid samples.
- The positive control, negative control, HBcAg Reagent, Anti-HBc·Peroxidase Conjugate Solution and specimens should be regarded as potential hazards to health. They should be used and discarded according to the user's laboratory safety procedures.
- Although all human sourced material are tested non-reactive for Anti-HCV

and Anti-HIV, and inactivated at 56°C for one hour, the reagent shall be handled as potential infectious material.

- Return any unused microplate strips to the plate pouch with desiccant.
- Prior to beginning the assay procedure, bring all reagents and required number of strips to room temperature (20-30°C).
- The reagent preparation method might be different from lot to lot, so please check the lot and follow the instructions given in this manual.
- Avoid air bubbles in the wells as this could result in lower binding efficiency and higher CV% of duplicate reading.
- Briefly spin down the all vials before use.
- If crystals are observed in the 20X Wash Buffer, warm to 37°C until the crystals are completely dissolved.
- Minimize lag time between wash steps to ensure the plate does not become completely dry during the assay.
- Take care not to contaminate the TMB Substrate. Do not expose the TMB solution to glass, foil or metal. If the solution is blue before use, DO NOT USE IT.
- Do NOT return leftover TMB Substrate to bottle. Do NOT contaminate the unused TMB Substrate. If the solution is blue before use, DO NOT USE IT.
- Change pipette tips between the addition of different reagent or samples.
- Taping the well strips together with lab tape can be done as an extra precaution to avoid plate strips from coming loose during the procedure.
- Run both standards and samples in at least duplicates (triplicate is recommended).

SAMPLE COLLECTION & STORAGE INFORMATION

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

<u>Serum</u>: Collect blood in a tube with no anticoagulant. Allow the blood to clot at room temperature for 30 minutes. Centrifuge at 2500 x g for 20 minutes at 4°C.

<u>**Plasma:**</u> Collect blood with EDTA, heparin or citrate and centrifuge at 2000 x g for 10 minutes at 4°C.

Note:

- Whole blood specimen should be separated as soon as possible in order to avoid hemolysis.
- Specimens must be stored at 4°C and avoided heat inactivation to minimize deterioration.
- For long term storage, specimens should be frozen below -20°C. Storage in self-defrosting freezers is not recommended.
- > Avoid multiple freeze-thaw procedures.
- The specimen must not contain any AZIDE compounds which can inhibit the peroxidase activity.
- Incompletely coagulated serum samples and microbial-contaminated specimens should not be used.

REAGENT PREPARATION

1X Wash Buffer: Dilute 20X Wash Buffer into distilled water to yield 1X
Wash Buffer. (E.g., add 25 mL of 20X Wash Buffer into 475 mL of distilled water to a final volume of 500 mL)

ASSAY PROCEDURE

All materials should be equilibrated to room temperature (RT, 20-30°C) before use. Standards and samples should be assayed in duplicates.

1. Make 1:100 dilution of each specimen. Prepare for each specimen a tube for dilution. Add 500 μ L of Specimen Diluent Buffer and 5 μ L of each specimen to each tube respectively and shake to mix.

Note: Do not dilute the controls.

- Reserve 1 wells for blanks. Add 100 μL of Negative Control to each three wells, 100 μL of Positive Control to each two wells, and 100 μL of Specimen Diluent Buffer to each other of the reaction wells for test sample.
- Add 5 μL of each diluted specimen to each well containing Specimen Diluent Buffer, respectively.
- 4. Gently tap the plate. Cover the plate and incubate for **1 hour** at **37°C**.
- Remove plate cover. Aspirate each well and wash, repeating the process 7 times for a total 8 washes. Wash by filling each well with **1X Wash Buffer** (350 μL) using a squirt bottle, manifold dispenser. Complete removal of liquid at each is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating, decanting or blotting against clean paper towels.

- 6. Add **50 μL** of **HBcAg Reagent** to each reaction well except the blank well.
- 7. Add **50 μL** of **Anti-HBc·Peroxidase Solution** to each well.
- 8. Gently tap the plate. Cover the plate and incubate for 1 hour at 37°C.
- 9. Aspirate each well and wash as step 5.
- Add 50 μL of TMB Substrate A first, then add 50 μL of TMB Substrate B to each well including the blank well. Mix well gently.
- 11. Cover the plate and incubate for **30 minutes** at **room temperature** in the dark.
- 12. Immediately Add $100 \,\mu$ L of Stop Solution to each well including the blank well. The color of the solution should change from blue to yellow.
- 13. Read the absorbance at **O.D. 450 nm**. (with **620-690 nm** reference wavelength)

Note: absorbance value of blank well must be less than 0.100.

CALCULATION OF RESULTS

1. Calculation of the NCx (Mean Absorbance of Negative Control).

Example:

Sample No. Absorbance

1	0.068
2	0.072
3	0.070

NCx = (0.068 + 0.072 + 0.070) / 3 = 0.070 (NCx must be ≤ 0.1 , otherwise,

the test run is invalid).

2. Calculation of the PCx (Mean Absorbance of Positive Control).

Example:	
Sample No.	Absorbance
1	1.612
2	1.613

PCx = (1.612 + 1.613) / 2 = 1.6125 (PCx must be ≥ 0.4 , otherwise, the test run is invalid).

3. Calculation of the P – N value

Example:

P - N = PCx - NCx = 1.6125 - 0.070 = 1.5425 (P-N value must be ≥ 0.3 ,

otherwise, the test run is invalid).

4. Calculation of the Cutoff Value

Cutoff Value = NCx + (PCx / 4)

Example: Cutoff Value =0.070 + (1.6125 / 4) = 0.473

5. Calculation of the Retest Range

Retest Range = Cutoff Value \pm 10%

Example: Cutoff Value = 0.473,

Retest Range = (0.473 - 0.047) to (0.473 + 0.047) = 0.426 to 0.520

INTERPRETATION OF RESULT

- Specimen with signal / Cut-off ration ≤ 0.9 are considered non-reactive for Anti-HBc IgM. Specimen with signal / Cut-off ration ≥ 1.1 are considered reactive for Anti-HBc IgM.
- If the signal / Cut-off ratio is within Retest Range, the test must be repeated in duplicate and interpreted as above. If both results are non-reactive the final result is non-reactive, if both results are reactive the final result is reactive. Any other combination is an indeterminate result. Testing of follow up samples and other hepatitis B serological markers should be taken into account in case of an indeterminate result.

QUALITY ASSURANCE

Sensitivity

2.5 PEI U/mL