

Ras Activation ELISA Kit

Ras Activation ELISA Kit is an Enzyme Immunoassay kit for the quantification Ras Activation cell sample.

Catalog number: ARG83384

Package: 96 wells

For research use only. Not for use in diagnostic procedures.

TABLE OF CONTENTS

SECTION	Page
INTRODUCTION	3
PRINCIPLE OF THE ASSAY	3
MATERIALS PROVIDED & STORAGE INFORMATION	4
MATERIALS REQUIRED BUT NOT PROVIDED	4
TECHNICAL NOTES AND PRECAUTIONS	5
REAGENT PREPARATION	6
SAMPLE COLLECTION & STORAGE INFORMATION	6
ASSAY PROCEDURE	7
CALCULATION OF RESULTS	義書籤。
EXAMPLE OF TYPICAL STANDARD CURVE	義書籤。
QUALITY ASSURANCE錯誤! 尚未定	義書籤。

MANUFACTURED BY:

Arigo Biolaboratories Corporation

Address: No. 22, Ln. 227, Gongyuan Rd., Hsinchu City 300, Taiwan

Phone: +886 (3) 562 1738

Fax: +886 (3) 561 3008

Email: info@arigobio.com

INTRODUCTION

The Ras superfamily, derived from "Rat sarcoma virus", is a protein superfamily of small GTPases. Members of the superfamily are divided into families and subfamilies based on their structure, sequence and function. The five main families are Ras, Rho, Ran, Rab and Arf GTPases. The Ras family itself is further divided into 6 subfamilies: Ras, Ral, Rap, Rheb, Rad and Rit. Miro is a recent contributor to the superfamily. Each subfamily shares the common core G domain, which provides essential GTPase and nucleotide exchange activity. The surrounding sequence helps determine the functional specificity of the small GTPase, for example the 'Insert Loop', common to the Rho subfamily, specifically contributes to binding to effector proteins such as WASP. In general, the Ras family is responsible for cell proliferation: Rho for cell morphology, Ran for nuclear transport, and Rab and Arf for vesicle transport.

PRINCIPLE OF THE ASSAY

This assay employs the quantitative sandwich enzyme immunoassay technique. An antibody specific for RAS has been pre-coated onto a microtiter plate. Samples are pipetted into the wells and any RAS present is bound by the immobilized antibody. After washing away any unbound substances, a biotinconjugated antibody specific for RAS is added to each well and incubate. Following a washing to remove unbound substances, streptavidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After washing away any unbound antibody-enzyme reagent, a substrate solution (TMB) is added to the wells and color develops in proportion to the amount of RAS bound in the initial step. The color development is

stopped by the addition of acid and the intensity of the color is measured at a wavelength of 450nm \pm 2nm.

MATERIALS PROVIDED & STORAGE INFORMATION

Upon received, store Raf-1 RBD at \leq -80°C, store 1000X Ras Activation Antibody and control at \leq -20°C, store other component at 2-8°C at all times.

Component	Quantity	Storage information
Antibody Coated Microplate	8 x 12 strips	4°C.
Raf-1 RBD	40 μL	≤ -80°C
1000X pan-Ras Antibody	1 vial (20 μL)	≤ -20°C
1000X HRP-Streptavidin Solution	1 vial (20 μL)	4°C
Diluent Buffer	50 mL	4°C
10X Wash Buffer	100 mL	4°C
5X Lysis Buffer	30 mL	4°C
TMB Substrate	12 mL	4°C (Protect from light)
Control	2 X 50 μL	≤ -20°C
Stop Solution	12 mL	4°C

MATERIALS REQUIRED BUT NOT PROVIDED

- Microplate reader capable of reading at O.D. 450 nm
- Centrifuge and centrifuge tube
- Deionized or Distilled water
- Pipettes, pipette tips and Multichannel micropipette reservoir
- Automated microplate washer (optional)

TECHNICAL NOTES AND PRECAUTIONS

- Wear protective gloves, clothing, eye, and face protection especially while handling blood or body fluid samples.
- Upon received, store Raf-1 RBD at ≤ -80°C, store 1000X Ras Activation
 Antibody and control at ≤ -20°C, store other component at 2-8°C at all
 times. Use the kit before expiration date.
- Prior to beginning the assay procedure, bring all reagents and required number of strips to room temperature.
- Briefly spin down the all vials before use.
- If crystals are observed in the 10X Wash Buffer, warm to 37°C until the crystals are completely dissolved.
- Minimize lag time between wash steps to ensure the plate does not become completely dry during the assay.
- Ensure complete reconstitution and dilution of reagents prior to use.
- Take care not to contaminate the TMB Substrate. Do not expose the TMB solution to glass, foil or metal. Do NOT return leftover TMB Substrate to bottle. Do NOT contaminate the unused TMB Substrate. If the solution is blue before use, DO NOT USE IT.
- Change pipette tips between the addition of different reagent or samples.
- Taping the well strips together with lab tape can be done as an extra
 precaution to avoid plate strips from coming loose during the procedure.

REAGENT PREPARATION

- 1X Wash Buffer: Dilute 10X Wash Buffer into distilled water to yield 1X
 Wash Buffer.
- **1X Lysis Buffer**: Dilute 5X Lysis Buffer into distilled water to yield **1X** Lysis Buffer. Add protease inhibitors just before use.
- 1X pan-Ras Antibody: 10 minutes before use, dilute 1000X pan-Ras Antibody into Diluent Buffer to yield 1X pan-Ras Antibody.
- 1X HRP-Streptavidin Solution: 10 minutes before use, dilute 2500X HRP-Streptavidin Solution into Diluent Buffer to yield 1X HRP-Streptavidin Solution. Keep diluted HRP-Streptavidin Solution in dark before use.

SAMPLE COLLECTION & STORAGE INFORMATION

The sample collection and storage conditions listed below are intended as general guidelines. Sample stability has not been evaluated.

Suspension Cells: Centrifuge $3-6 \times 10^6$ cells at $1200 \times g$ for 2 minutes and discard supernatant. Wash cell pellet twice with ice-cold PBS, centrifuge, and discard the supernatant. Resuspend cell pellet in 0.5-1 mL of cold 1X Lysis Buffer. Lyse cells with repeated pipetting on ice. Centrifuge at $12000 \times g$ for 10 minutes and collect the cell lysate supernatant.

Adherent Cells: Wash $1-5 \times 10^6$ cells twice with ice-cold PBS per 100 mm dish, and add 0.5 - 1 mL of ice-cold 1X Lysis Buffer for 10 - 20 min. Harvest cells with a cell scraper. Centrifuge at $14000 \times g$ for 10 minutes at $4^{\circ}C$ and collect the cell lysate supernatant.

ASSAY PROCEDURE

All materials should be equilibrated to room temperature (RT, 20-25°C) before use. All assay should be assayed in duplicates.

Control Loading

- 1. Add 0.5 mL same samples to two microcentrifuge tubes.
- 2. Add 10 µL 0.5 M EDTA to each sample.
- 3. Add 5 μ L of positive control and 5 μ L of 100X negative control to the other tube.
- 4. Mix and incubate for 30 minutes at 30°C.
- 5. Add 33 μ L of 1 M MgCl2 to each tube for stopping the loading.
- 6. Mix and place tubes on ice.

Ras Activation ELISA step

- 1. Determine the number of wells to be used, and dilute the Raf-1 RBD 1:500 in Assay Diluent. Add **100 µL** of **1X** Raf-1 RBD into each wells the plate.
- 2. Cover the plate and incubate for **1 hour** at **RT**.
- 3. Aspirate each well and wash, repeating the process 2 time for a **total 3** washes. Wash by filling each well with **1X Wash Buffer (250 μL)** using a squirt bottle, manifold dispenser, or autowasher. Complete removal of liquid at each is essential to good performance. After the last wash, remove any remaining Wash Buffer by aspirating, decanting or blotting against clean paper towels.
- 4. Add **50 μL** of **lysate sample**, **control** and **blank** to well.
- 5. Add **50 \muL** of **Diluent Buffer** to each well (100 μ L total volume).
- 6. Cover the plate and incubate for **1 hour** at **RT**.

Ras Activation ELISA Kit ARG83384

- 7. Aspirate each well and wash plate as step 3, but for 5 5 times wash.
- 8. Add $100 \mu L$ of 1X pan-Ras Antibody to each well.
- 9. Cover the plate and incubate for **1 hour** at **RT**.
- 10. Aspirate each well and wash plate as step 7.
- 11. Add **100 μL** of **1X HRP-Streptavidin Solution** to each well.
- 12. Cover the plate and incubate for **1 hour** at **room temperature** in the dark.
- 13. Aspirate each well and wash plate as step 7.
- 14. Add 100 μL of TMB Substrate in each well.
- 15. Incubate for **5-20 mins** at **room temperature** in the dark.
- 16. Add **100 μL** of **Stop Solution** to each well to stop the reaction.
- 17. Read the absorbance with a plate reader at **O.D. 450 nm.** It is recommended reading the absorbance within 10 minutes after adding the stop solution.