

ARG11046 anti-Hepatocellular Carcinoma antibody [CHALV1]

Package: 100 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [CHALV1] recognizes Hepatocellular Carcinoma
Tested Reactivity	Hu
Tested Application	ICC/IF
Specificity	There were no detectable positive reactions with other human tumor cell lines and solid tumour specimens. Showed negative reaction to all normal human adult tissues and fetal tissues tested. Using indirect immunofluorescence assay, this mAb showed positive reactions with 8 human hepato-cellular carcinoma cell lines and 12 human liver cancer surgical specimens tested.
Host	Mouse
Clonality	Monoclonal
Clone	CHALV1
Isotype	IgG1, kappa
Target Name	Hepatocellular Carcinoma
Species	Human
Immunogen	Human liver carcinoma cells.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ICC/IF	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS (pH 7.0)
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.