

ARG11047
anti-CA50 antibody [93]Package: 100 µg
Store at: -20°C

Summary


Product Description	Mouse Monoclonal antibody [93] recognizes CA50
Tested Reactivity	Hu
Tested Application	ELISA, ICC/IF
Specificity	React with human CA50.
Host	Mouse
Clonality	Monoclonal
Clone	93
Isotype	IgG3, kappa
Target Name	CA50
Species	Human
Immunogen	Human CA50 protein.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	0.2 ng/ml
	ICC/IF	10 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS (pH 7.2)
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG11047 anti-CA50 antibody [93] ICC/IF image

Immunofluorescence: Methanol-fixed HeLa cells stained with ARG11047 anti-CA50 antibody [93] at 10 µg/ml dilution.
