

## ARG11161 anti-Rabies virus Glycoprotein antibody [7E3]

Package: 200 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [7E3] recognizes Rabies virus Glycoprotein
Tested Reactivity	Virus
Tested Application	EIA, ELISA, ICC/IF, WB
Specificity	MAbs 7E3 are specific to Rabies virus glycoprotein.
Host	Mouse
Clonality	Monoclonal
Clone	7E3
Isotype	IgG2a
Target Name	Rabies virus Glycoprotein
Immunogen	purified rabies virus
Conjugation	Un-conjugated

### Application Instructions

Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	--

### Properties

Form	Liquid
Purification	Protein A affinity purified.
Buffer	PBS (pH 7.4) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Concentration	1.0-2.0 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Microbiology and Infectious Disease antibody
---------------	--