

ARG11168 anti-ACAN / Aggrecan antibody

Package: 50 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes ACAN / Aggrecan
Tested Reactivity	Ms, Rat
Predict Reactivity	Hu, Bov, Chk, Dog, Hrs, Pig, Rb
Tested Application	IHC-P
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	ACAN / Aggrecan
Species	Human
Conjugation	Un-conjugated
Alternate Names	AGC1; SEDK; AGCAN; CSPG1; MSK16; CSPGCP

Application Instructions

Application table	Application	Dilution
	IHC-P	1:100 - 1:500

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties


Form	Liquid
Purification	Purification with Protein A.
Buffer	Aqueous buffered solution, 0.09% Sodium azide, 50% Glycerol and 1% BSA.
Preservative	0.09% Sodium azide
Stabilizer	50% Glycerol and 1% BSA
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	ACAN
-------------	------

Gene Full Name aggrecan
Calculated Mw 261 kDa

Images



ARG11168 anti-ACAN / Aggrecan antibody IHC-P image

Immunohistochemistry: Rat synovial stained with ARG11168 anti-ACAN / Aggrecan antibody.