

## ARG20720 Rat IgG2b Isotype Control antibody [KLH/G2b-1-2] (Cyanine 5)

Package: 50 µg  
Store at: 4°C

### Summary

Product Description	Cyanine 5-conjugated Rat Monoclonal antibody [KLH/G2b-1-2] as a negative control antibody for Rat IgG2b
Tested Application	ELISA, FACS, FLISA, IHC-Fr, IHC-P
Specificity	KLH
Host	Rat
Clonality	Monoclonal
Clone	KLH/G2b-1-2
Isotype	IgG2b, kappa
Target Name	Rat IgG2b
Conjugation	Cyanine 5

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 0.1 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.