

ARG21570 Goat anti-Mouse IgG2c antibody (HRP), pre-adsorbed

Package: 500 µl
Store at: 4°C

Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Mouse IgG2c
Tested Reactivity	Ms
Tested Application	Depletion, ELISA, ELISPOT, FACS, FLISA, IHC-P, WB
Specificity	Reacts with the heavy chain of C57BL/6 mouse IgG2c
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG2c
Species	Mouse
Conjugation	HRP

Application Instructions

Pre Adsorbed Mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA; Human immunoglobulins and pooled sera.

Application table


Application	Dilution
Depletion	Assay-dependent
ELISA	1:4000 - 1:8000
ELISPOT	Assay-dependent
FACS	Assay-dependent
FLISA	Assay-dependent
IHC-P	Assay-dependent
WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	50% PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Images



ARG21570 Goat anti-Mouse IgG2c antibody (HRP) (pre-adsorbed) ELISA image

ELISA: The plate was coated with purified Mouse IgG1, IgG2a, IgG2b, IgG2c, IgG3, IgM, and IgA. Immunoglobulins were detected with serially diluted ARG21570 Goat anti-Mouse IgG2c antibody (HRP) (pre-adsorbed).