

ARG21759 Rabbit anti-Rat IgG (H+L) antibody, F(ab')₂ fragment (AP)

Package: 500 µl
Store at: 4°C

Summary

Product Description	AP-conjugated F(ab') ₂ fragment of Rabbit Polyclonal antibody recognizes Rat IgG (H+L)
Tested Reactivity	Rat
Tested Application	ELISA, FLISA, ICC/IF, IHC-Fr, IHC-P, WB
Specificity	Reacts with the heavy and light chains of rat IgG and the light chains of rat IgM
Host	Rabbit
Clonality	Polyclonal
Isotype	F(ab') ₂ IgG
Target Name	IgG (H+L)
Species	Rat
Conjugation	AP

Application Instructions

Application table	Application	Dilution
	ELISA	1:2000 - 1:4000
	FLISA	Assay-dependent
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	50 mM Tris (pH 8.0), 0.1% Sodium azide, 1 mM MgCl ₂ and 50% Glycerol
Preservative	0.1% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.