

ARG21902 Goat anti-Mouse IgG3 antibody (Biotin), pre-adsorbed

Package: 500 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Goat Polyclonal antibody recognizes Mouse IgG3
Tested Reactivity	Ms
Tested Application	ELISA, ELISPOT, FACS, FLISA, ICC/IF, IHC-Fr, IHC-P, WB
Specificity	Reacts with the heavy chain of mouse IgG3
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG3
Species	Mouse
Conjugation	Biotin

Application Instructions

Pre Adsorbed Mouse IgG1, IgG2a, IgG2b, IgM, and IgA; Human immunoglobulins and pooled sera.

Application table

Application	Dilution
ELISA	1:5000 - 1:20000
ELISPOT	Assay-dependent
FACS	< 1 µg/10 ⁶ cells
FLISA	Assay-dependent
ICC/IF	Assay-dependent
IHC-Fr	Assay-dependent
IHC-P	Assay-dependent
WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml


Storage instruction

Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.

Images



ARG21902 Goat anti-Mouse IgG3 antibody (Biotin) (pre-adsorbed) ELISA image

ELISA: The plate was coated with purified Mouse IgG1, IgG2a, IgG2b, IgG3, IgM, and IgA. Immunoglobulins were detected with serially diluted ARG21902 Goat anti-Mouse IgG3 antibody (Biotin) (pre-adsorbed) followed by [ARG23912](#) Streptavidin (HRP).