

**ARG22313**  
**anti-Influenza A nucleoprotein antibody [FluA-NP 4F1] (HRP)**Package: 200 µl  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Mouse Monoclonal antibody [FluA-NP 4F1] recognizes Influenza A nucleoprotein
Tested Reactivity	Influenza A
Tested Application	ELISA, ICC/IF, WB
Host	Mouse
Clonality	Monoclonal
Clone	FluA-NP 4F1
Isotype	IgG1
Target Name	Influenza A nucleoprotein
Species	Virus
Immunogen	Recombinant influenza virus type A nucleoprotein
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	ELISA	1:2000 - 1:8000
	ICC/IF	Assay-dependent
	WB	1:2000 - 1:8000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	50% PBS (pH 7.4) and 50% Glycerol.
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.