

ARG22685 anti-Chlamydia trachomatis LPS antibody [1312/236]

Package: 100 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [1312/236] recognizes Chlamydia trachomatis LPS This antibody recognizes C. trachomiatitis lipopolysaccharide. In a simple ELISA this antibody is reactive with 15 serovars of C. trachomatis. The epitope has been shown by Western blotting to be present in LPS.
Tested Reactivity	Bacteria
Tested Application	ELISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	1312/236
Isotype	IgG1
Target Name	Chlamydia trachomatis LPS
Species	Bacteria
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purification with Protein A.
Buffer	PBS and 0.09% Sodium azide.
Preservative	0.09% Sodium azide
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.