

ARG22688
anti-HIV1 p24 antibodyPackage: 1 ml
Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes HIV1 p24 This antibody recognizes HIV1 viral capsid p24 protein encoded by the gag gene along with p6, p7 and p17 which provide structural elements of the virus. Goat anti Human immunodeficiency virus 1 p24 antibody does not cross-react with Bovine sera, Human sera or Human T or B cells.
Tested Reactivity	HIV
Tested Application	ICC/IF, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	HIV1 p24
Species	HIV
Immunogen	Purified native p24, from strain IIIB.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ICC/IF	Assay-dependent
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Unpurified.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	4 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.