

**ARG22699**  
**Goat anti-Cat IgG (Fc) antibody**Package: 500 µg  
Store at: -20°C

### Summary

Product Description	Goat Polyclonal antibody recognizes Cat IgG (Fc) This antibody recognizes feline IgG (Fc) and shows no cross-reactivity with other feline immunoglobulin classes in immunoelectrophoresis.
Tested Reactivity	Cat, Rb
Tested Application	ELISA
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (Fc)
Species	Cat
Immunogen	Purified feline IgG (Fc) fragment.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	1:100 - 1:1000

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.09% Sodium azide.
Preservative	0.09% Sodium azide
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.