

**ARG22816**  
**anti-V5 tag antibody [SV5-Pk1] (HRP)**Package: 50 µg  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Mouse Monoclonal antibody [SV5-Pk1] recognizes V5 tag V5-Tag antibody, clone SV5-Pk1 recognizes the sequence, IPNPLLGL, present on the P/V proteins of the paramyxovirus, SV5. It detects recombinant proteins, including transmembrane and secreted proteins tagged with this sequence.
Tested Reactivity	Other
Tested Application	ELISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	SV5-Pk1
Isotype	IgG2a
Target Name	V5 tag
Species	Others
Immunogen	Paramyxovirus Simian-Virus 5 (SV5)
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purification with Protein A.
Buffer	PBS, 0.01% Thiomersal and HRP Stabiliser
Preservative	0.01% Thiomersal
Stabilizer	HRP Stabiliser
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.