

## ARG22817 anti-6X His tag antibody [AD1.1.10] (Biotin)

Package: 50 µg  
Store at: 4°C

### Summary

Product Description	Biotin-conjugated Mouse Monoclonal antibody [AD1.1.10] recognizes 6X His tag
Tested Reactivity	Other
Tested Application	ELISA, WB
Host	Mouse
Clonality	Monoclonal
Clone	AD1.1.10
Isotype	IgG1
Target Name	6X His tag
Species	Others
Immunogen	PAX6 transcription factor linked to histidine tag.
Conjugation	Biotin

### Application Instructions


Application table	Application	Dilution
	ELISA	Assay-dependent
	WB	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS and 0.09% Sodium azide
Preservative	0.09% Sodium azide
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

## Images



ARG22817 anti-6X His tag antibody [AD1.1.10] (Biotin) WB image

Western blot: HIS tagged protein stained with ARG22817 anti-6X His tag antibody [AD1.1.10] (Biotin).