

## ARG23729 Rabbit anti-Hamster IgG (H+L) antibody (HRP), pre-adsorbed

Package: 500 µl  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Rabbit Polyclonal antibody recognizes Hamster IgG (H+L)
Tested Reactivity	Hm
Species Does Not React With	Ms
Tested Application	ELISA, FACS, FLISA
Specificity	Reacts with the heavy and light chains of Syrian and Armenian hamster IgG.
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Hamster
Immunogen	Hamster IgG (H+L)
Conjugation	HRP

### Application Instructions

Pre Adsorbed Mouse immunoglobulins and pooled sera.

#### Application table

Application	Dilution
ELISA	1:2000 - 1:4000
FACS	Assay-dependent
FLISA	Assay-dependent

Application Note \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.