

ARG23730 Rabbit anti-Hamster IgG (H+L) antibody (HRP), pre-adsorbed

Package: 500 µl
Store at: 4°C

Summary

Product Description	HRP-conjugated Rabbit Polyclonal antibody recognizes Hamster IgG (H+L)
Tested Reactivity	Hm
Species Does Not React With	Ms, Rat
Tested Application	ELISA, FACS, FLISA, IHC-Fr, WB
Specificity	Reacts with the heavy and light chains of Armenian and Syrian hamster IgG as demonstrated by ELISA or FLISA. May also react with other species and with the light chains of other hamster immunoglobulins. Minimal cross reactivity with Mouse and Rat serum proteins.
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Hamster
Immunogen	Hamster IgG (H+L)
Conjugation	HRP

Application Instructions

Pre Adsorbed Mouse and Rat serum proteins. May react with immunoglobulins from other species and the light chains of other hamster immunoglobulins.

Application table	Application	Dilution
	ELISA	1:4000 - 1:8000
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.4) and 50% Glycerol
Stabilizer	50% Glycerol

Storage instruction

Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.