

ARG23813 Goat anti-Mouse IgG2b antibody (Spectral Red), pre-adsorbed

Package: 125 µg
Store at: 4°C

Summary

Product Description	Spectral Red-conjugated Goat Polyclonal antibody recognizes Mouse IgG2b
Tested Reactivity	Ms
Tested Application	ELISA, ELISPOT, FACS, FLISA, ICC/IF, IHC-Fr, IHC-P, WB
Specificity	Reacts with the heavy chain of mouse IgG2b
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG2b
Species	Mouse
Conjugation	Spectral Red

Application Instructions

Pre Adsorbed Mouse IgG1, IgG2a, IgG3, IgM, and IgA; Human immunoglobulins and pooled sera.

Application table

Application	Dilution
ELISA	Assay-dependent
ELISPOT	Assay-dependent
FACS	< 0.1 µg/10 ⁶ cells
FLISA	Assay-dependent
ICC/IF	Assay-dependent
IHC-Fr	Assay-dependent
IHC-P	Assay-dependent
WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose

Concentration	0.25 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.