

ARG23822 Goat anti-Rat IgG (H+L) antibody, pre-adsorbed

Package: 250 µg
Store at: -20°C

Summary

Product Description	Goat Polyclonal antibody recognizes Rat IgG (H+L)
Tested Reactivity	Rat
Tested Application	ELISA, FACS
Specificity	Reacts with the heavy and light chains of rat IgG and the light chains of rat IgM
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Rat
Conjugation	Un-conjugated

Application Instructions

Pre Adsorbed Mouse and Human immunoglobulins and pooled sera.

Application table

Application	Dilution
ELISA	Assay-dependent
FACS	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	BBS (pH 8.2).
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.