

ARG23829 Mouse anti-Rat Kappa Light Chain antibody [K4F5]

Package: 250 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [K4F5] recognizes Rat Kappa Light Chain
Tested Reactivity	Rat
Tested Application	ELISA, FACS, FLISA, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	K4F5
Isotype	IgG2a, kappa
Target Name	Kappa Light Chain
Species	Rat
Immunogen	Rat Kappa hybridoma
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	Assay-dependent
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	BBS (pH 8.2).
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.