

## ARG23833 Mouse anti-Rat Kappa Light Chain antibody [K4F5] (Biotin)

Package: 250 µg  
Store at: 4°C

### Summary

Product Description	Biotin-conjugated Mouse Monoclonal antibody [K4F5] recognizes Rat Kappa Light Chain
Tested Reactivity	Rat
Tested Application	ELISA, FACS, FLISA, IHC-Fr, WB
Host	Mouse
Clonality	Monoclonal
Clone	K4F5
Isotype	IgG2a, kappa
Target Name	Kappa Light Chain
Species	Rat
Immunogen	Rat Kappa hybridoma
Conjugation	Biotin

### Application Instructions

Application table	Application	Dilution
	ELISA	1:5000 - 1:10000
	FACS	< 1 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-Fr	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.