

## ARG23850 Mouse anti-Rabbit IgG antibody [2A9] (Spectral Red)

Package: 50 µg  
Store at: 4°C

### Summary

Product Description	Spectral Red-conjugated Mouse Monoclonal antibody [2A9] recognizes Rabbit IgG
Tested Reactivity	Rb
Tested Application	ELISA, FACS, FLISA, IHC-P, WB
Host	Mouse
Clonality	Monoclonal
Clone	2A9
Isotype	IgG1
Target Name	IgG
Species	Rabbit
Conjugation	Spectral Red

### Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	FACS	< 0.3 µg/10 <sup>6</sup> cells
	FLISA	Assay-dependent
	IHC-P	Assay-dependent
	WB	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.