

## ARG23919 Goat anti-Monkey IgG (H+L) antibody (HRP)

Package: 500 µg  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Monkey IgG (H+L)
Tested Reactivity	Hu, AGMK, Bb, Chimp, Mk, R. Mk
Tested Application	ELISA, ID
Specificity	This antibody has been shown to react specifically with monkey (African green monkey, Cynomolgus monkey, Baboon, Chimpanzee and Rhesus Monkey) IgG heavy chains and with light chains common to other monkey immunoglobulins by immunoelectrophoresis and ELISA. Reaction to human was also detected. The antibody may cross react with IgG from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Goat
Immunogen	Purified monkey IgG.
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	ELISA	1:10000 - 1:100000
	ID	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS, 0.05% Proclin 300 and 0.2% BSA.
Preservative	0.05% Proclin 300
Stabilizer	0.2% BSA
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.