

## ARG23926 anti-Aspergillus antibody

Package: 500 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes Aspergillus
Tested Reactivity	Fungi
Tested Application	ELISA
Specificity	This antibody recognizes soluble proteins from common Aspergillus species. The antibody has not been cross absorbed and may react with related microorganisms.
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Aspergillus
Species	Fungi
Immunogen	Soluble extract from A. fumigatus, A. flavus, A. niger and A. terreus.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ELISA	1:100 - 1:1500
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purified.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	4 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.