

ARG24040
anti-T-2 Toxin antibody [BCAB1360]Package: 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [BCAB1360] recognizes T-2 Toxin
Tested Reactivity	Other
Tested Application	ELISA
Specificity	Cross Reactivity: T-2 Toxin: 100% HT-2 Toxin: 117%
Host	Mouse
Clonality	Monoclonal
Clone	BCAB1360
Isotype	IgG1
Target Name	T-2 Toxin
Species	Others
Immunogen	HT-2 Toxin-Human serum albumin immune conjugate.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purification with Protein G.
Purity	> 95% (by electrophoresis)
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.