

ARG24041
anti-Streptomycin antibody [BCAB1331]Package: 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [BCAB1331] recognizes Streptomycin
Tested Reactivity	Other
Tested Application	ELISA
Specificity	The antibody reacts with Streptomycin. The antibody does not cross react with Human serum albumin.
Host	Mouse
Clonality	Monoclonal
Clone	BCAB1331
Isotype	IgG1
Target Name	Streptomycin
Species	Others
Immunogen	Streptomycin-Human serum albumin immune conjugate.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purification with Protein G.
Purity	> 95% (by electrophoresis)
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.