

ARG24046 anti-Aflatoxin B2 antibody [BCAB1312]

Package: 50 µg
Store at: -20°C

Summary

Product Description	Mouse Monoclonal antibody [BCAB1312] recognizes Aflatoxin B2
Tested Reactivity	Other
Tested Application	ELISA
Specificity	Cross Reactivity: Aflatoxin B1: 100% Aflatoxin B2: 29% Aflatoxin G2: 6% Aflatoxin M1: 5% Aflatoxin G1: 1%
Host	Mouse
Clonality	Monoclonal
Clone	BCAB1312
Isotype	IgG1
Target Name	Aflatoxin B2
Species	Others
Immunogen	Aflatoxin B2.
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purification with Protein G.
Purity	> 95% (by electrophoresis)
Buffer	Neutralized citrate buffer (pH 7.4) and 0.05% Sodium azide.
Preservative	0.05% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.