

## ARG24084 Goat anti-Chicken IgA antibody (HRP)

Package: 500 µg  
Store at: 4°C

### Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Chicken IgA
Tested Reactivity	Chk
Tested Application	ELISA
Specificity	This antibody shows no cross-reactivity with other chicken immunoglobulin classes in immunoelectrophoresis. This antibody may react with IgA from other species.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgA
Species	Chicken
Immunogen	Purified chicken IgA.
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	ELISA	1:10000 - 1:100000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS, 0.05% Proclin 300 and 0.2% BSA.
Preservative	0.05% Proclin 300
Stabilizer	0.2% BSA
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.