

ARG24179 anti-Asialofetuin antibody

Package: 1 ml
Store at: -20°C

Summary

Product Description	Sheep Polyclonal recognizes Asialofetuin
Tested Reactivity	Bov
Tested Application	ELISA
Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Target Name	Asialofetuin
Immunogen	Native
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS and 0.09% Sodium Azide
Preservative	0.09% Sodium Azide
Concentration	5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.