

ARG30049 SyK / Zap70 Antibody Duo

Package: 1 pair
Store at: -20°C

Component

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG62632	anti-Syk antibody [4D10.1]	Mouse mAb	Hu	FACS, IHC-P, IP, WB	250 µl
ARG65322	anti-ZAP70 antibody	Goat pAb	Hu	IP, WB	25 µg

Summary


Target Name	SyK / Zap70
Alternate Names	SyK / Zap70 antibody; Spleen tyrosine kinase / Zeta chain of T cell receptor associated protein kinase 70 Antibody; Syk antibody; ZAP70 antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.


Bioinformation

Gene Full Name	Spleen tyrosine kinase (SyK) / Zeta chain of T cell receptor associated protein kinase 70 (ZAP70) Antibody Duo
Research Area	Cell Biology and Cellular Response antibody; Controls and Markers antibody; Immune System antibody; Signaling Transduction antibody



ARG65322 anti-ZAP70 antibody WB image

Western Blot: Jurkat lysate (35 µg protein in RIPA buffer) stained with ARG65322 anti-Zap70 antibody at 0.3 µg/ml dilution.



ARG65322 anti-ZAP70 antibody WB image

Western blot: 30 µg of Jurkat and Molt4 cell lysates stained with ARG65322 anti-ZAP70 antibody at 1:500 dilution.