

ARG30291
serum Amyloid A ELISA Antibody DuoPackage: 1 pair
Store at: -20°C**Component**

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10639	anti-serum Amyloid A antibody [SAA19]	Mouse mAb	Cat	ELISA, WB	50 µg
ARG10640	anti-serum Amyloid A antibody [VSA34]	Mouse mAb	Cat	ELISA, WB	50 µg

Summary

Target Name	serum Amyloid A
Alternate Names	serum Amyloid A ELISA antibody; serum Amyloid A antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.