

ARG30348 Exosome Antibody Panel

Package: 1 kit
Store at: -20°C

Component


Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG57952	anti-CD63 antibody	Rabbit pAb	Hu	ICC/IF, IHC-P, WB	20 µl
ARG40187	anti-TSG101 antibody	Rabbit pAb	Hu, Ms, Rat	FACS, ICC/IF, IHC-P, WB	20 µl
ARG57583	anti-GM130 antibody	Rabbit pAb	Hu, Ms, Rat, Cow, Dog, Mk	ICC/IF, IHC-P, IP, WB	20 µl

Properties

Storage instruction For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.


Note For laboratory research only, not for drug, diagnostic or other use.

Images



ARG57952 anti-CD63 antibody IHC-P image


Immunohistochemistry: Paraffin-embedded Mouse testis stained with ARG57952 anti-CD63 antibody at 1:100 dilution.



A375

ARG57952 anti-CD63 antibody WB image


Western blot: 25 µg of A375 cell lysate stained with ARG57952 anti-CD63 antibody at 1:1000 dilution.



Jurkat


ARG40187 anti-TSG101 antibody WB image

Western blot: Jurkat cell lysate stained with ARG40187 anti-TSG101 antibody.



ARG57583 anti-GM130 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Human cervix carcinoma stained with ARG57583 anti-GM130 antibody.



MCF7

ARG57583 anti-GM130 antibody WB image

Western blot: MCF7 cell lysate stained with ARG57583 anti-GM130 antibody.