

## ARG40048 anti-TIPE2 antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes TIPE2
Tested Reactivity	Hu, Ms, Rat
Tested Application	ICC/IF, IP, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	TIPE2
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 1-184 of Human TIPE2 (NP_078851.2).
Conjugation	Un-conjugated
Alternate Names	Inflammation factor protein 20; TNF alpha-induced protein 8-like protein 2; Tumor necrosis factor alpha-induced protein 8-like protein 2; TNFAIP8-like protein 2; TIPE2

### Application Instructions

Application table	Application	Dilution
	ICC/IF	1:50 - 1:200
	IP	Assay-dependent
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Mouse thymus	
Observed Size	20 kDa	

### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

**Note**

For laboratory research only, not for drug, diagnostic or other use.

## Bioinformation

---

Gene Symbol	TNFAIP8L2
Gene Full Name	tumor necrosis factor, alpha-induced protein 8-like 2
Function	Acts as a negative regulator of innate and adaptive immunity by maintaining immune homeostasis. Negative regulator of Toll-like receptor and T-cell receptor function. Prevents hyperresponsiveness of the immune system and maintains immune homeostasis. Inhibits JUN/AP1 and NF-kappa-B activation. Promotes Fas-induced apoptosis (By similarity). [UniProt]
Calculated Mw	21 kDa

## Images

---



ARG40048 anti-TIPE2 antibody WB image

Western blot: 25 µg of Mouse thymus lysate stained with ARG40048 anti-TIPE2 antibody at 1:1000 dilution.