

ARG41308 anti-HA tag antibody (Biotin)

Package: 100 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Rabbit Polyclonal antibody recognizes HA tag
Tested Reactivity	Other
Tested Application	ELISA, IE, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	HA tag
Species	Others
Immunogen	KLH-conjugated HA peptide.
Conjugation	Biotin

Application Instructions

Application table	Application	Dilution
	ELISA	0.2 - 0.05 µg/ml
	IE	
	WB	0.2 - 1 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.2), 0.1% Sodium azide and 0.2% BSA.
Preservative	0.1% Sodium azide
Stabilizer	0.2% BSA
Concentration	1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.