

## ARG41953 anti-6X His tag antibody (HRP)

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	HRP-conjugated Mouse Monoclonal antibody recognizes 6X His tag
Tested Reactivity	Other
Tested Application	WB
Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Target Name	6X His tag
Species	Others
Immunogen	A synthetic peptide of His tag
Conjugation	HRP

### Application Instructions

Application table	Application	Dilution
	WB	1:2000 - 1:20000


**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Affinity Purified
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Controls and Markers antibody; Tag Internal Control antibody
---------------	--



ARG41953 anti-6X His tag antibody (HRP)

Western blot: 20 ng of His-tagged protein stained with ARG41953 anti-6X His tag antibody (HRP) at 1:20,000 dilution.