

ARG42167
anti-HOXC9 antibodyPackage: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes HOXC9
Tested Reactivity	Zfish
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	HOXC9
Species	Zebrafish
Immunogen	KLH-conjugated synthetic peptide between aa. 234-260 of Zebrafish HOXC9.
Conjugation	Un-conjugated
Alternate Names	HOX3; Homeobox protein Hox-3B; Homeobox protein Hox-C9; HOX3B; Hox-C9; Homeobox protein Hox-C9a; bz39; BZ-39; hoxc9

Application Instructions

Application table	Application	Dilution
	WB	1:1000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	Zebrafish muscle	
Observed Size	~ 30 kDa	


Properties

Form	Liquid
Purification	Purification with Protein A and immunogen peptide.
Buffer	PBS and 0.09% (W/V) Sodium azide.
Preservative	0.09% (W/V) Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	hoxc9a
Gene Full Name	homeobox C9a
Background	This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXC genes located in a cluster on chromosome 12. [provided by RefSeq, Jul 2008]
Function	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. [UniProt]
Calculated Mw	30 kDa
Cellular Localization	Nucleus. [UniProt]

Images



ARG42167 anti-HOXC9 antibody WB image

Western blot: 20 µg of Zebrafish muscle lysate stained with ARG42167 anti-HOXC9 antibody at 1:1000 dilution.