

ARG42327
anti-CD68 antibody [Y1/82A] (APC)

Package: 50 tests

Store at: 4°C

Summary

Product Description	APC-conjugated Mouse Monoclonal antibody [Y1/82A] recognizes CD68
Tested Reactivity	Hu
Tested Application	FACS
Specificity	The mouse monoclonal antibody Y1/82A recognizes CD68 (LAMP4), a 110 kDa glycoprotein expressed mainly in cytoplasmic granules of monocytes/macrophages, granulocytes, and dendritic cells.
Host	Mouse
Clonality	Monoclonal
Clone	Y1/82A
Isotype	IgG2a
Target Name	CD68
Species	Human
Immunogen	Lysosomal contents of lung macrophages.
Conjugation	APC
Alternate Names	Macrosialin; CD antigen CD68; LAMP4; Gp110; GP110; SCARD1

Application Instructions

Application table	Application	Dilution
	FACS	10 µl / 100 µl of whole blood or 10 ⁶ cells

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.


Properties

Form	Liquid
Purification	Purified
Buffer	PBS and 15 mM Sodium azide.
Preservative	15 mM Sodium azide
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	CD68
Gene Full Name	CD68 molecule
Background	This gene encodes a 110-kD transmembrane glycoprotein that is highly expressed by human monocytes and tissue macrophages. It is a member of the lysosomal/endosomal-associated membrane glycoprotein (LAMP) family. The protein primarily localizes to lysosomes and endosomes with a smaller fraction circulating to the cell surface. It is a type I integral membrane protein with a heavily glycosylated extracellular domain and binds to tissue- and organ-specific lectins or selectins. The protein is also a member of the scavenger receptor family. Scavenger receptors typically function to clear cellular debris, promote phagocytosis, and mediate the recruitment and activation of macrophages. Alternative splicing results in multiple transcripts encoding different isoforms. [provided by RefSeq, Jul 2008]
Function	Could play a role in phagocytic activities of tissue macrophages, both in intracellular lysosomal metabolism and extracellular cell-cell and cell-pathogen interactions. Binds to tissue- and organ-specific lectins or selectins, allowing homing of macrophage subsets to particular sites. Rapid recirculation of CD68 from endosomes and lysosomes to the plasma membrane may allow macrophages to crawl over selectin-bearing substrates or other cells. [UniProt]
Highlight	Related products: CD68 antibodies ; CD68 Duos / Panels ; Anti-Mouse IgG secondary antibodies ; Related news: Exploring Antiviral Immune Response Anti-SerpinB9 therapy, a new strategy for cancer therapy RIP1 activation and pathogenesis of NASH
Calculated Mw	37 kDa
PTM	N- and O-glycosylated. [UniProt]
Cellular Localization	Isoform Short: Cell membrane; Single-pass type I membrane protein. Isoform Long: Endosome membrane; Single-pass type I membrane protein. Lysosome membrane; Single-pass type I membrane protein. [UniProt]


Images



ARG42327 anti-CD68 antibody [Y1/82A] (APC) FACS image

Flow Cytometry: Human peripheral whole blood stained with ARG42327 anti-CD68 antibody [Y1/82A] (APC) at 10 μ l / 100 μ l of peripheral whole blood.

ARG42327 anti-CD68 antibody [Y1/82A] (APC) FACS image



Flow Cytometry: Separation of Human monocytes (red-filled) from Human CD68 negative lymphocytes (black-dashed). Human peripheral whole blood stained with ARG42327 anti-CD68 antibody [Y1/82A] (APC) at 10 μ l / 100 μ l of peripheral whole blood.