

ARG43900 anti-FOXO3A antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes FOXO3A
Tested Reactivity	Hu
Tested Application	ELISA, WB
Host	Rabbit
Clonality	Polyclonal
Target Name	FOXO3A
Species	Human
Immunogen	Human FOXO3A recombinant protein
Conjugation	Un-conjugated
Alternate Names	FOXO3; Forkhead Box O3; AF6q21; FKHL1; FOXO3A; FOXO2; Forkhead In Rhabdomyosarcoma-Like 1; Forkhead Box Protein O3; Forkhead, Drosophila, Homolog Of, In Rhabdomyosarcoma-Like 1; Forkhead Homolog (Rhabdomyosarcoma) Like 1; Forkhead Box O3A; AF6q21 Protein; FKHL1P2

Application Instructions

Application table	Application	Dilution
	ELISA	0.1-0.5 µg/ml
	WB	0.25-0.5 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties


Form	Liquid
Buffer	0.9% NaCl, 0.2% Na ₂ HPO ₄ and 4% Trehalose.
Stabilizer	4% Trehalose
Concentration	0.5 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation

Gene Symbol	FOXO3
-------------	-------

Gene Full Name	Forkhead Box O3
Background	This gene belongs to the forkhead family of transcription factors which are characterized by a distinct forkhead domain. This gene likely functions as a trigger for apoptosis through expression of genes necessary for cell death. Translocation of this gene with the MLL gene is associated with secondary acute leukemia. Alternatively spliced transcript variants encoding the same protein have been observed.
Function	Transcriptional activator that recognizes and binds to the DNA sequence 5'-[AG]TAAA[TC]A-3' and regulates different processes, such as apoptosis and autophagy
Calculated Mw	71 kDa
PTM	Acetylation, Methylation, Phosphoprotein, Ubl conjugation
Cellular Localization	Cytoplasm, Membrane, Mitochondrion, Mitochondrion outer membrane, Nucleus

Images



ARG43900 anti-FOXO3A antibody WB image

Western blot: PC-3 and 293T stained with ARG43900 anti-FOXO3A antibody at 0.5 µg/mL dilution.