

**ARG44324**  
anti-Rattin antibodyPackage: 50 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes Rattin
Tested Reactivity	Rat
Tested Application	WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Rattin
Species	Rat
Immunogen	Synthetic peptide
Conjugation	Un-conjugated
Alternate Names	HNr; Humanin-like protein; Rattin

### Application Instructions

Application table	Application	Dilution
	WB	1:500-1:1000

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Purification	Antigen Affinity Purified
Buffer	PBS with 0.02% Sodium azide
Preservative	0.02% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

### Bioinformation

Gene Symbol	HNr
Gene Full Name	Humanin-like protein
Function	Plays a role as a neuroprotective factor.