

ARG44845 anti-Histone H3.1 + Histone H3.3 antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rat Monoclonal antibody recognizes Histone H3.1 (P68431) and Histone H3.3 (P84243).
Tested Reactivity	Hu
Tested Application	IHC-P, IP, WB
Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Target Name	Histone H3.1 / Histone H3.3
Conjugation	Un-conjugated

Application Instructions

Application table	Application	Dilution
	IHC-P	1:500
	IP	1:100
	WB	1:1000

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Protein A purification
Buffer	PBS with 0.09% sodium azide
Preservative	0.09% sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation

Gene Symbol	H3C1; H3-3A
Gene Full Name	H3 Clustered Histone 1; H3.3 Histone A