

ARG45007 anti-IgA antibody [RM220]

Package: 50 µg
Store at: -20°C

Summary

Product Description	ARG45007 Rabbit Monoclonal antibody [RM220] recognizes IgA.
Tested Reactivity	Ms
Tested Application	ELISA, WB
Specificity	This antibody specifically reacts to mouse IgA. No cross reactivity with mouse IgG1, IgG2a, IgG2b, IgG2c, IgG3, IgM, IgE, or human IgA.
Host	Rabbit
Clonality	Monoclonal
Clone	RM220
Isotype	IgG
Target Name	IgA
Species	Mouse
Immunogen	Mouse IgA
Target Ig	IgA
Conjugation	Un-conjugated
Alternate Names	Immunoglobulin A

Application Instructions

Application table	Application	Dilution
	ELISA	0.005 - 0.2 µg/mL
	WB	0.5 - 2 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Buffer	PBS with 50% Glycerol, 1% BSA and 0.09% sodium azide
Preservative	0.09% sodium azide
Stabilizer	50% Glycerol, 1% BSA and 0.09%
Concentration	1 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Note

For laboratory research only, not for drug, diagnostic or other use.