

ARG54490 anti-6X His tag antibody (HRP)

Package: 50 µg
Store at: 4°C

Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes 6X His tag
Tested Reactivity	Other
Tested Application	ELISA, ICC/IF, IHC, WB
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	6X His tag
Species	Others
Immunogen	6-His (HHHHHH) coupled to KLH. Antibody was isolated by affinity chromatography using the peptide immobilized on solid support then conjugated to horseradish peroxidase (HRP). NOTE: This antibody does not detect all N-terminal 6-His tagged proteins.
Conjugation	HRP

Application Instructions

Application Note	ELISA: 1:10,000 - 1:100,000 Western blot: Colorimetric: 1:1,000 - 1:10,000 Chemiluminescent: 1:1,000 - 1:30,000 Immunohistochemistry: not recommended Immunocytochemistry: 1:200 - 1:500 * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Properties

Form	Liquid
Buffer	PBS, 0.2% BSA and 0.1% Pro-Clean 400
Preservative	0.1% Pro-Clean 400
Stabilizer	0.2% BSA
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Controls and Markers antibody; Tag Internal Control antibody
---------------	--------------------------------------------------------------