

## ARG54491 anti-6X His tag antibody

Package: 100 µg  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes 6X His tag
Tested Reactivity	Other
Tested Application	ELISA, ICC/IF, IHC, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	6X His tag
Species	Others
Immunogen	6-His (HHHHHH) coupled to KLH.
Conjugation	Un-conjugated

### Application Instructions

Application Note	ELISA: Coating: 2 - 10 µg/ml Primary: 0.03 - 1 µg/ml Western blot: Colorimetric: 0.1 - 1 µg/ml Chemiluminescent: 0.03 - 1 µg/ml IHC, ICC/IF: 0.5 - 5 µg/ml * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.
------------------	---

### Properties

Form	Liquid
Buffer	PBS (pH 7.2) and 0.1% Sodium azide
Preservative	0.1% Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Research Area	Controls and Markers antibody; Tag Internal Control antibody
---------------	--