

ARG55095 anti-HA tag antibody - CHIP-Grade

Package: 100 µl
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes HA Tag - CHIP-Grade
Tested Reactivity	Other
Tested Application	ChIP, ChIP-seq, ICC/IF, IP, R-IP A, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	HA tag
Species	Others
Immunogen	KLH-conjugated encoding HA tag (CYPYDVPDYAYPYDVPDYA)
Conjugation	Un-conjugated

Application Instructions


Application table	Application	Dilution
	ChIP	1:100
	ChIP-seq	Assay-dependent
	ICC/IF	1:10 - 1:50
	IP	Assay-dependent
	R-IP A	Assay-dependent
	WB	1:1000

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS and 0.09% (W/V) Sodium azide
Preservative	0.09% (W/V) Sodium azide
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG55095 anti-HA tag antibody - CHIP-Grade WB image

Western blot: 10 μ g of 293T and 293T + myc-HA-mCherry cell lysates stained with ARG55095 anti-HA tag antibody - CHIP-Grade at 1:1000 dilution.