

## ARG56177 anti-EBV Early Antigen antibody [1108-1]

Package: 50 µg  
Store at: -20°C

### Summary

Product Description	Mouse Monoclonal antibody [1108-1] recognizes EBV Early Antigen
Tested Reactivity	EBV
Tested Application	ICC/IF
Host	Mouse
Clonality	Monoclonal
Clone	1108-1
Isotype	IgG1, kappa
Target Name	EBV Early Antigen
Species	EBV
Immunogen	Affinity purified early antigen polypeptides from induced Raji cells precipitated with African Burkitt's lymphoma sera.
Conjugation	Un-conjugated

### Application Instructions

Application table	Application	Dilution
	ICC/IF	0.5 - 1 µg/ml
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Purification	Purification with Protein G.
Buffer	PBS (pH 7.4), 0.05% Sodium azide and 0.1 mg/ml BSA.
Preservative	0.05% Sodium azide
Stabilizer	0.1 mg/ml BSA
Concentration	0.2 mg/ml
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.