

ARG56575 anti-Helicobacter pylori Vacuolating cytotoxin autotransporter antibody

Package: 50 µg
Store at: -20°C

Summary

Product Description	Rabbit Polyclonal antibody recognizes Helicobacter pylori Vacuolating cytotoxin autotransporter
Tested Reactivity	H. pylori
Tested Application	ELISA
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	Helicobacter pylori Vacuolating cytotoxin autotransporter
Species	H. pylori
Immunogen	Recombinant H. pylori Vacuolating cytotoxin autotransporter protein.
Conjugation	Un-conjugated
Alternate Names	vacA

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

Properties

Form	Liquid
Purification	Purified
Buffer	PBS (pH 7.4), 0.03% Proclin-300 and 50% Glycerol.
Preservative	0.03% Proclin-300
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight

Related products:

[Helicobacter pylori antibodies](#); [Helicobacter pylori ELISA Kits](#); [Helicobacter pylori Duos / Panels](#); [Anti-Rabbit IgG secondary antibodies](#);

Related news:

[Tools for studying H. pylori diseases](#)

Research Area

Microbiology and Infectious Disease antibody